ARCTIC SUN[®]
TEMPERATURE MANAGEMENT SYSTEM

ARCTICGEL[™] Small Universal Pad

(REF 318-01)

INSTRUCTIONS FOR USE

ARCTICGEL™ Small Universal Pad - Instructions for Use**Intended Use**

The ARCTIC SUN® ARCTICGEL™ Pads are intended for use with the ARCTIC SUN® Temperature Management System, to provide energy (heat) transfer between the patient and the temperature-controlled water circulating through the ARCTICGEL™ Pads in order to provide targeted temperature management.

Indications for Use

The ARCTIC SUN® Temperature Management System is a thermal regulating system, indicated for monitoring and controlling patient temperature in adult and pediatric patients of all ages.

Contraindications

- There are no known contraindications for the use of a thermoregulatory system.
- Do not place ARCTICGEL™ Pads on skin that has signs of ulcerations, burns, hives or rash.
- Do not place ARCTICGEL™ Pads on immature (non-keratinized) skin or premature babies.
- While there are no known allergies to hydrogel materials, caution should be exercised with any patient with a history of skin allergies or sensitivities.

Warning

Do not place ARCTICGEL™ Pads over transdermal medication patches as warming can increase drug delivery, resulting in possible harm to the patient.

Cautions

- Federal law restricts this device to sale by or on the order of a physician.
- The ARCTICGEL™ Small Universal Pads are only for use with the ARCTIC SUN® Temperature Management System.
- This product is to be used by or under the supervision of trained, qualified medical personnel.
- The clinician is responsible for determining the appropriateness of use of this device and the user-settable parameters, including water temperature, for each patient.
- Due to underlying medical or physiological conditions, some patients are more susceptible to skin damage from pressure and heat or cold. Patients at risk include those with poor tissue perfusion or poor skin integrity due to edema, diabetes, peripheral vascular disease, poor nutritional status or steroid or high dose vasopressor therapy. If accessible, examine the patient's skin under the ARCTICGEL™ Pads often; especially those patients at higher risk of skin injury.
- Skin injury may occur as a cumulative result of pressure, time and temperature. Possible skin injuries include bruising, tearing, skin ulcerations, blistering, and necrosis. Do not place bean bags or other firm positioning devices under the ARCTICGEL™ Pads. Do not place any positioning devices under the pad manifolds or patient lines.
- If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury.
- Do not allow urine, antibacterial solutions or other agents to pool underneath the ARCTICGEL™ Pads. Urine and antibacterial agents can absorb into the pad hydrogel and cause chemical injury and loss of pad adhesion. Replace pads immediately if these fluids come into contact with the hydrogel.
- Do not place ARCTICGEL™ Pads directly over an electrosurgical grounding pad. The combination of heat sources may result in skin burns.
- If needed, place defibrillation pads between the ARCTICGEL™ Pads and the patient's skin.
- The ARCTICGEL™ Pads are non-sterile for single patient use only. Do not reprocess or sterilize. If used in a sterile environment, pads should be placed according to the physician's request, either prior to the sterile preparation or sterile draping.
- Use pads immediately after opening. Do not store pads in opened pouch.
- The ARCTICGEL™ Pads should not be punctured with sharp objects. Punctures will result in air entering the fluid pathway and may reduce performance.
- The ARCTICGEL™ Pads are for single patient use. The water content of the hydrogel affects the pad's adhesion to the

skin and conductivity, and therefore, the efficiency of controlling patient temperature. Periodically check that pads remain moist and adherent. Replace pads when the hydrogel no longer uniformly adheres to the skin. Replacing pads at least every 5 days is recommended.

- Carefully and slowly remove ARCTICGEL™ Pads from the patient's skin at the completion of use. Aggressive removal of the pad from the patient's skin may result in skin tears.
- Do not allow circulating water to contaminate the sterile field when lines are disconnected.
- Discard used ARCTICGEL™ Pads in accordance with hospital procedures for medical waste.

Directions

- The ARCTICGEL™ Small Universal Pads are only for use with the ARCTIC SUN® Temperature Management System. See Operators Manual for detailed instructions on system use.
- Select the proper number of pads for the patient size and clinical indication. However, the rate of temperature change and potentially the final achievable temperature is affected by pad surface area, patient size, pad placement and water temperature range. Best system performance will be achieved by using the maximum number of pads.
- The following sizing chart is provided as guidance. Modify the number of pads and locations as necessary to meet specific clinical needs.

Patient Weight	Number of Pads
2.5 - 5 kg (5.5 - 10 lbs)	1
5 - 10 kg (10 - 22 lbs)	2
10 - 16 kg (22 - 35 lbs)	3
16 - 30 kg (35 - 66 lbs)	Use ARCTICGEL™ Pad Kit 317-02 (XXS)

- Place the pads on healthy, clean skin only. Locate pad edges away from articulating areas of the body to avoid irritation. Place pads to allow for full respiratory excursion.
- The pad surface must be contacting the skin for optimal energy transfer efficiency. The Small Universal Pads are provided with a cloth liner over the hydrogel. The pads may be used with the cloth liner in place or removed to expose the hydrogel adhesive. If the cloth liner is used, secure the pads with the Velcro straps to ensure good contact with the skin. The pads may be removed and reapplied if necessary.
- Due to the small patient size and the potential for rapid patient temperature change
 - It is recommended to use the following settings:
Water Temperature High Limit: ≤40°C (104°F)
Water Temperature Low Limit: ≥10°C (50°F)
Control Strategy: 2
 - It is recommended to use the Patient Temperature High and Patient Temperature Low alert settings.
- Place a core patient temperature probe and connect to the ARCTIC SUN® Temperature Management System patient temperature 1 connector for continuous patient temperature feedback. A rectal or esophageal temperature probe is recommended.
- Attach the pad's line connectors to the fluid delivery line manifolds. Begin treating the patient.
- If the pads fail to prime or a significant continuous air leak is observed in the pad return line, check connections, then if needed replace the leaking pad.
- Once all the pads are primed, assure the steady state flow rate displayed on the control panel is appropriate for the number of pads connected.

Number Pads	1	2	3
Minimum Flow Rate l/m	0.9	1.7	2.4

- When finished, purge water from pads. Slowly remove pads from the patient and discard.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Symbols Glossary

The FDA required English symbols glossary for Becton Dickinson (BD) is located at <https://www.bd.com/en-us/symbols-glossary>.

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are trademarks and/or registered
trademarks of C. R. Bard, Inc.
© 2022 C. R. Bard, Inc. All rights reserved.
Medivance, Inc.
Becton Dickinson Ireland Ltd.